





GN32 Powder Free Latex Exam Gloves

White

Latex

Powder free $\sqrt{}$

Textured 🗹

Gloves tested to comply with Medical Devices Directives 93/42/EEC and European Standards EN455 parts 1, 2 and 3, EN388 and EN374. Designed for medical procedures, offering a high quality and dependable latex glove.

 $\overline{\mathbf{A}}$

- Visually inspected for holes and minor defects.
- Random samples are filled with 1000ml of water and examined for leaks.
- Tested to AQL 1.5

M

- Manufactured from natural rubber latex.
- Typical force at break : 9 Newtons.
- Typical elongation : 700%
- Typical tensile strength : 21 Mpa.
- Available from size XS to XL



Dimensions	Extra Small	Small	Medium	Large	Extra Large
Palm Width (mm)	≤ 80 mm	80 ±10 mm	95 ±10 mm	110 ±10 mm	≥ 110 mm
Length (mm)	Min 240 mm	Min 240 mm	Min 240 mm	Min 240 mm	Min 240 mm
Thickness (mm) palm	0.11 ± 0.02	0.11 ± 0.02	0.11 ± 0.02	0.11 ± 0.02	0.11 ± 0.02
Thickness (mm) finger	0.14 ± 0.02	0.14 ± 0.02	0.14 ± 0.02	0.14 ± 0.02	0.14 ± 0.02
Thickness (mm) cuff	0.08 ± 0.02	0.08 ± 0.02	0.08 ± 0.02	0.08 ± 0.02	0.08 ± 0.02

Benefits		
Low protein	Minimising the risk of skin irritation and allergic reactions	
Low in residual	With no detectable thiurams, reduces the risk of contact	
chemicals	dermatitis associated with these chemicals	
Textured	To improve grip	
Powder free and chlorinated	To reduce dust contamination and latex protein content	
International testing	Tested to Lowry Method, Endotoxin and Chemical Residual	
Beaded cuff	For added strength and easier donning	
Available in 5 sizes	Providing wearer with a range of comfortable sizes	

Characteristics	Unaged	Aged	
Force at break	9 Newtons	6 Newtons	
Elongation at break	700%	650%	

Recommended Use

Medical, Laboratory, Dental, Beauty, Veterinary, Pharmaceutical

HPC Healthline Limited Colwood House 158 Garth Road Morden, Surrey England, UK, SM4 4LZ


T 0044 (0) 20 8335 3636

F 0044 (0) 20 8335 0300

E sales@hpchealthline.co.uk

W www.hpchealthline.co.uk

GN32 Powder Free Latex Exam Gloves





Storage instructions Y Keep in a cool dry place Keep away from direct sunlight or fluorescent light Store at temperature no less than 5°C, no more than 25°C Shelf life : 3 years

Packing

100 gloves per dispenser

10 dispensers per case 72 cases per pallet Morden, Surrey England, UK , SM4 4LZ T 0044 (0) 20 8335 3636

F 0044 (0) 20 8335 0300

HPC Healthline Limited Colwood House

158 Garth Road

E sales@hpchealthline.co.uk

W www.hpchealthline.co.uk



Fine Touch